

# Anti-phospho-FAK (Tyr576) Rabbit antibody



Catalog Number: 167690

## Product name

Anti-phospho-FAK (Tyr576) Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

phospho-FAK (Tyr576) Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Tyr576 of human FAK

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications


WB, FC

## Dilutions

WB: 1/2000-1/10000

FC: 1/20

## Validations



Western blot detection of phospho-FAK (Tyr576) in A549,U251,C2C12 cell lysates using phospho-FAK (Tyr576) Rabbit mAb(1:1000 diluted).Predicted band size:119kDa.Observed band size:125kDa.