

Anti-TRIM24 Rabbit antibody



Catalog Number: 167692

Product name

Anti-TRIM24 Rabbit antibody

Specificity

Human, Mouse

Antibody description

TRIM24 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human TRIM24

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/1000

IHC: 1/50

ICC/IF: 1/20

Validations



Western blot detection of TRIM24 in HeLa,A549 cell lysates using TRIM24 Rabbit mAb(1:1000 diluted).Predicted band size:117kDa.Observed band size:117kDa.