

Anti-RPA70 Rabbit antibody



Catalog Number: 167802

Product name

Anti-RPA70 Rabbit antibody

Specificity

Human

Antibody description

RPA70 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human RPA70

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/2000-1/10000

FC: 1/20

IP: 1/20

Validations



Western blot detection of RPA70 in A549,HL-60 cell lysates using RPA70 Rabbit mAb(1:1000 diluted).Predicted band size:68kDa.Observed band size:70kDa.