

Anti-Galectin 1 Rabbit antibody



Catalog Number: 167869

Product name

Anti-Galectin 1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Galectin 1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Galectin 1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of Galectin 1 in CHO-K1, C6 cell lysates using Galectin 1 Rabbit mAb(1:1000 diluted). Predicted band size: 15kDa. Observed band size: 15kDa.