

Anti-IKB epsilon Rabbit antibody



Catalog Number: 167895

Product name

Rabbit IgG

Anti-IKB epsilon Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF, IP

Human

Dilutions

Antibody description

WB: 1/1000

IKB epsilon Rabbit monoclonal antibody

IHC: 1/50

Preparation

ICC/IF: 1/50

Antigen: A synthetic peptide of human IKB epsilon

IP: 1/20

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Validations

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Western blot detection of IKB epsilon in HeLa, Ramos using IKB epsilon Rabbit mAb(1:1000 diluted)

