

Anti-GRK2 Rabbit antibody



Catalog Number: 167945

Product name

Anti-GRK2 Rabbit antibody

Specificity

Human

Antibody description

GRK2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human GRK2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB: 1/1000

IP: 1/20

Validations



Western blot detection of GRK2 in HeLa, A549, HL-60, U2OS cell lysates using GRK2 Rabbit mAb(1:1000 diluted). Predicted band size: 80kDa. Observed band size: 80kDa.