

Anti-TATA binding protein TBP Rabbit antibody



Catalog Number: 167989

Product name

Anti-TATA binding protein TBP Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TATA binding protein TBP Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human TATA binding protein TBP

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/500

ICC/IF: 1/50

FC: 1/50

IP: 1/20

Validations



Western blot detection of TATA binding protein TBP in HeLa,A549,HL-60,U2OS cell lysates using TATA binding protein TBP Rabbit mAb(1:1000 diluted).Predicted band size:38kDa.Observed band size:38kDa.