

# Anti-VAV2 Rabbit antibody



Catalog Number: 168030

## Product name

Anti-VAV2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

VAV2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human VAV2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/5000

IHC: 1/20

ICC/IF: 1/50

IP: 1/20

## Validations



Western blot detection of VAV2 in HeLa, MCF-7, HepG2, U2OS, U251 using VAV2 Rabbit mAb(1:1000 diluted)