

Anti-Osteopontin Rabbit antibody



Catalog Number: 168038

Product name

Anti-Osteopontin Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Osteopontin Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Osteopontin

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC

Validations



Western blot detection of Osteopontin in HeLa, CHO-K1, C6, Jurkat cell lysates using Osteopontin Rabbit mAb(1:1000 diluted). Predicted band size: 35kDa. Observed band size: 35kDa.