

# Anti-DTYMK Rabbit antibody



Catalog Number: 168053

## Product name

Anti-DTYMK Rabbit antibody

## Specificity

Human

## Antibody description

DTYMK Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human DTYMK

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/50

FC: 1/20

IP: 1/20

## Validations



Western blot detection of DTYMK in HeLa,Ramos cell lysates using DTYMK Rabbit mAb(1:1000 diluted).Predicted band size:24kDa.Observed band size:24kDa.