

# Anti-CDX2 (4F12) Mouse antibody



Catalog Number: 168114

## Product name

Anti-CDX2 (4F12) Mouse antibody

## Specificity

Human, Monkey

## Antibody description

CDX2 (4F12) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CDX2 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, IP, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/100

## Validations




Immunocytochemistry stain of HeLa using CDX2



mouse mAb (1:100).

Western blot detection of CDX2 in HeLa, COS7 and K562 cell lysates using CDX2 mouse mAb (1:1000 diluted). Predicted band size: 34KDa. Observed



band size: 34KDa.

Immunoprecipitation analysis of HeLa cell lysate using CDX2 mouse mAb.