

Anti-TXNRD2 (3F2) Mouse antibody



Catalog Number: 168243

Product name

Anti-TXNRD2 (3F2) Mouse antibody

Specificity

Human

Antibody description

TXNRD2 (3F2) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human TXNRD2 protein fragments expressed in E.coli.

Formulation

PBS(pH 7.4)containing with 0.02% NaN₃ 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of TXNRD2 in A431, U251, U2OS, Lncap and HeLa cell lysates and using TXNRD2 mouse mAb (1:1000 diluted). Predicted band size: 56KDa. Observed band size: 56KDa.