

# Anti-CKMT2 (3F4) Mouse antibody



Catalog Number: 168271

## Product name

Anti-CKMT2 (3F4) Mouse antibody

## Specificity

Rat

## Antibody description

CKMT2 (3F4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CKMT2 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of CKMT2 in Rat Brain lysates using CKMT2 mouse mAb (1:1000 diluted). Predicted band size: 47KDa. Observed band size: 47KDa.