

# Anti-CD133 (8F2) Mouse antibody



Catalog Number: 168404

## Product name

Anti-CD133 (8F2) Mouse antibody

## Specificity

Human

## Antibody description

CD133 (8F2) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human Prominin-1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB

## Dilutions

WB: 1/3000

## Validations



Western blot detection of CD133 in CaCo2 cell lysate using CD133 mouse mAb (1:3000 diluted). Predicted band size: 133 kDa. Observed band size: 133 kDa