

Anti-Hsp60 (6C8) Mouse antibody



Catalog Number: 168407

Product name

Anti-Hsp60 (6C8) Mouse antibody

Specificity

Human, Mouse, Rat

Antibody description

Hsp60 (6C8) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Hsp60 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications

WB, ICC/IF, IP


Dilutions

WB: 1/1000


;ICC: 1/100

;IP: 1/100


Validations



Western blot detection of Hsp60 in HeLa, Raji, Jurkat, C6, C2C12 and 3T3 cell lysates using Hsp60 mouse mAb (1:1000 diluted). Predicted band size: 60kDa. Observed band size: 60kDa.



Immunocytochemistry staining of HeLa cells fixed with -20u2103 Methanol and using Hsp60 mouse



mAb (dilution 1:100).

Immunoprecipitation analysis of HeLa cell lysates using Hsp60 mouse mAb.