

# Anti-NEDD8 (5B8) Mouse antibody



Catalog Number: 168420

## Product name

Anti-NEDD8 (5B8) Mouse antibody

## Specificity

Human

## Antibody description

NEDD8 (5B8) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human NEDD8 expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

## Validations



Western blot detection of NEDD8 in C6 and HeLa cell lysates using NEDD8 mouse mAb (dilution 1:500). Predicted band size: 9kDa. Observed band size: 90kDa.