

Anti-GRK2 (6C1) Mouse antibody



Catalog Number: 168433

Product name

Anti-GRK2 (6C1) Mouse antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

GRK2 (6C1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human GRK2 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF

Dilutions


WB: 1/500 - 1/2000

IHC: 1/200 - 1/1000


ICC: 1/200 - 1/1000

ELISA: 1/10000


Validations



Western blot detection of GRK2 in C6, Jurkat, K562 and 3T3 cell lysates using GRK2 mouse mAb (1:1000 diluted). Predicted band size: 80KDa. Observed band size: 80KDa.



Immunohistochemical analysis of paraffin-embedded endometrial cancer tissues (left) and cervical cancer tissues (right) using GRK2 mouse mAb with DAB staining.



Immunofluorescence analysis of NIH/3T3 cells using GRK2 mouse mAb (green). Blue: DRAQ5 fluorescent DNA dye.