

Anti-GRAP2 (5C7) Mouse antibody



Catalog Number: 168444

Product name

Anti-GRAP2 (5C7) Mouse antibody

Ig Type

Mouse IgG2b

Specificity

Human

Applications

WB

Antibody description

GRAP2 (5C7) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Recombinant protein of human GRAP2.

Validations

Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal



Western blot detection of GRAP2 in GRAP2 positive cell (CEM, Molt-4, Jurkat) and negative cell (3T3, K562, HeLa) lysates using GRAP2 mouse mAb (1:1000 diluted). Predicted band size: 38KDa. Observed band size: 38KDa.