

Anti-ERCC1 (7F6) Mouse antibody



Catalog Number: 168445

Product name

Anti-ERCC1 (7F6) Mouse antibody

Specificity

Human

Antibody description

ERCC1 (7F6) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human ERCC1 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide 0.1mg/mlBSA and 50% glycerol.

Storage

Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB

Dilutions

1:1000

Validations



Western blot detection of ERCC1 in T47D,ZR75-1,CaCO2 and Molt-4 cell lysate using ERCC1 mouse mAb (1:1000 diluted).Predicted band size: 39KDa.Observed band size: 39KDa.