

# Anti-CDX2 (4B9) Mouse antibody



Catalog Number: 168446

## Product name

Anti-CDX2 (4B9) Mouse antibody

## Specificity

Human, Monkey, Mouse

## Antibody description

CDX2 (4B9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CDX2 protein fragments expressed in E.coli

## Formulation

Purified mouse monoclonal in PBS(pH 7.4)containing with 0.02% sodium azide 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


WB, IP, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/100

## Validations




Immunocytochemistry stain of HeLa using CDX2



mouse mAb (1:100).

Immunoprecipitation analysis of HeLa cell lysate using CDX2 mouse mAb.



Western blot detection of CDX2 in HeLa,Caco2,C2C12,SW480,COS7,HepG2 and 3T3 cell lysates using CDX2 mouse mAb (1:1000 diluted).Predicted band size:34KDa.Observed band size:34KDa.