

# Anti-CRK (3H7) Mouse antibody



Catalog Number: 168453

## Product name

Anti-CRK (3H7) Mouse antibody

## Specificity

Human

## Antibody description

CRK (3H7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human CrkII protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of CrkII in Lncap,A431,MCF7 and HT1080 cell lysates using CrkII mouse mAb (1:1000 diluted).Predicted band size: 34kDa.Observed band size: 34kDa.