

Anti-Skp1 (4E11) Mouse antibody



Catalog Number: 168539

Product name

Anti-Skp1 (4E11) Mouse antibody

Specificity

Human

Antibody description

Skp1 (4E11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human SKP1 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P, ICC/IF, FC

Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000

ICC: 1/200 - 1/1000

FC: 1/200 - 1/400

ELISA: 1/10000

Validations



Western blot detection of Skp1 in COS7, Raji, 3T3, PC12 and HeLa cell lysates using Skp1 mouse mAb (dilution 1:1000). Predicted band size: 19KDa. Observed band size: 19KDa.