

Anti-Ezrin Mouse antibody



Catalog Number: 168548

Product name

Anti-Ezrin Mouse antibody

WB: 1/500 - 1/2000

ICC: 1/200 - 1/1000

Specificity

Human

FC: 1/200 - 1/400

ELISA: 1/10000

Antibody description

Ezrin Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human EZR (AA: 292-464) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type


Mouse IgG1

Applications


WB, IHC-P, ICC/IF, FC

Dilutions

Validations



Western blot detection of Ezrin in Jurkat, A549, MCF7 and HeLa cell lysates using Ezrin mouse mAb (1:2000 diluted). Predicted band size: 69 kDa. Observed band size: 81 kDa.



Immunoprecipitation analysis of Jurkat cell lysates using Ezrin mouse mAb.