

Anti-PRPF4 Rabbit antibody



Catalog Number: 168634

Product name

Anti-PRPF4 Rabbit antibody

Specificity

Human

Antibody description

PRPF4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human PRPF4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/2000

ICC/IF: 1/200

FC: 1/20

IP: 1/20

Validations



Western blot detection of PRPF4 in K562,3T3,HeLa cell lysates using PRPF4 Rabbit mAb(1:1000 diluted).Predicted band size:59kDa.Observed band size:59kDa.