

Anti-FADD Rabbit antibody



Catalog Number: 168656

Product name

Anti-FADD Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse, Rat

Applications

WB, FC, IP

Antibody description

FADD Rabbit monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic peptide of human FADD

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of FADD in Rat Brain, 3T3 cell lysates using FADD Rabbit mAb (1:1000 diluted). Predicted band size: 23 kDa. Observed band size: 28 kDa.