

Anti-CDC7 Rabbit antibody



Catalog Number: 168733

Product name

Anti-CDC7 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CDC7 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human CDC7

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/100

IP: 1/20

Validations



Western blot detection of CDC7 in HeLa,CHO-K1,Jurkat cell lysates using CDC7 Rabbit mAb(1:500 diluted).Predicted band size:64kDa.Observed band size:64kDa.