

Anti-ABCG1 Rabbit antibody



Catalog Number: 168819

Product name

Anti-ABCG1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ABCG1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human ABCG1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

1:2000

1:20

Validations



Western blot detection of ABCG1 in HeLa cell lysates using ABCG1 Rabbit mAb(1:1000 diluted).Predicted band size:76kDa.Observed band size:100kDa.