

# Anti-GNAQ Rabbit antibody

Catalog Number: 168854

## Product name

Anti-GNAQ Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GNAQ Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human GNAQ

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

1:1000

## Validations



Western blot detection of GNAQ in Rat Brain, C6, 3T3 cell lysates using GNAQ Rabbit mAb(1:1000 diluted). Predicted band size: 42kDa. Observed band size: 42kDa.