

Anti-yars2 antibody



Catalog Number: 175176

Product name

Anti-yars2 antibody

Specificity

Zebrafish

Antibody description

Rabbit polyclonal antibody to yars2

Preparation

This antigen of this antibody was recombinant protein within zebrafish yars2 aa 271-494.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

IF, IHC-P

Dilutions

IF:1:50-1:200

IHC-P:1:50-1:200

Validations




Fig1:: IF staining of YARS2 in Zebrafish embryos (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.




Fig2:: Immunohistochemical analysis of paraffin-embedded Zebrafish embryos tissue section using anti-YARS2 antibody. Counter stained with hematoxylin.