

# Anti-Tnfrsf13b antibody



Catalog Number: 175729

## Product name

Anti-Tnfrsf13b antibody

## Specificity

Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to Tnfrsf13b

## Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from mouse tnfrsf13b 21-120/249

## Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

## Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

## Validations




Fig1: Sample:; Spleen (Mouse) Lysate at 40 ug; Lymph node (Mouse) Lysate at 40 ug; Primary: Anti- TNFRSF13B at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 32 kD; Observed band size: 30 kD