

Eribulin Monoclonal Antibody



Catalog Number: 176729

Product name

Eribulin Monoclonal Antibody

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

Eribulin conjugated with BSA. Protein G affinity purified

Formulation

PBS(pH7.4), lyophilized. **Reconstitution: Adding sterile water, prepare a stock solution of 2.0 mg/ml.**

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1, κ


Applications

ELISA

Dilutions


ELISA: 1:5,000-1:20,000

Validations



ELISA Binding Assay of Anti-Eribulin Antibody (MM461B) to different proteins.

ELISA analysis of Anti-Eribulin Monoclonal Antibody with Eribulin-BSA, Eribulin-OVA, BSA and OVA. Anti-Eribulin Monoclonal Antibody can bind to Eribulin-BSA or Eribulin-OVA but not carrier proteins.



ELISA Binding Assay of Anti-Eribulin Antibody (MM461B).

Immobilized Eribulin ADC at 2 μ g/mL (100 μ L/well) can bind anti-Eribulin mouse monoclonal antibody.