

Budesonide Monoclonal Antibody



Catalog Number: 176733

Product name

Budesonide Monoclonal Antibody

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

Budesonide conjugated with BSA. Protein G affinity purified

Formulation

PBS(pH7.4)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1, κ


Applications

ELISA

Dilutions


ELISA: 1:5,000-1:20,000

Validations



ELISA Binding Assay of Anti- Budesonide Antibody (MM817B) to different proteins.

ELISA analysis of Anti-Budesonide Monoclonal Antibody with Budesonide-BSA, Budesonide-OVA, BSA, OVA and Budesonide ADC. Anti-Budesonide Monoclonal Antibody can bind to Budesonide-BSA, Budesonide-OVA or Budesonide ADC but not carrier proteins.



ELISA Binding Assay of Anti- Budesonide Antibody (MM817B)

Immobilized Budesonide ADC at 2 μ g/ mL (100 μ L/well) can bind anti-Budesonide mouse monoclonal antibody.