

# Mouse CD123/IL3RA (Fc Tag) recombinant protein



Catalog Number: 500779

## General Information

### Gene Name Synonym

Interleukin-3 receptor class II alpha chain

### Protein Construction

A DNA sequence encoding the mouse Il3ra (NP\_032395.1) (Met1-Lys331) was expressed with the Fc region of human IgG1 at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Ser 17

## Molecular Mass

The recombinant mouse Il3ra consists 556 amino acids and predicts a molecular mass of 61.5 kDa.

## Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

## SDS-PAGE



Mouse CD123 / IL3RA Protein (Fc Tag)