

Human RELA / Transcription factor p65 (aa 1-306, GST Tag) recombinant protein



Catalog Number: 501219

General Information

Gene Name Synonym

Nuclear factor NF-kappa-B p65 subunit; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3

Protein Construction

A DNA sequence encoding the human RELA (Q04206-1) RHD domain (Met 1-Tyr 306) was fused with the GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human RELA/GST chimera

consists of 540 amino acids and has a predicted molecular mass of 62 kDa. It migrates as an approximately 58 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 0.15M NaCl, 20mM GST, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human RELA / Transcription factor p65 / NFkB p65 Protein (aa 1-306, GST Tag) SDS-PAGE