

# Human RAB31/RAB22B (His Tag) recombinant protein



Catalog Number: 501785

## General Information

### Gene Name Synonym

Ras-related protein Rab-22B

### Protein Construction

A DNA sequence encoding the human RAB31 (NP\_006859.2) (Met1-Cys195) was expressed with a polyhistidine tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

His

### Molecular Mass

The recombinant human RAB31 consists of 213

amino acids and predicts a molecular mass of 23.9 kDa. It migrates as an approximately 26 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 50 mM Tris, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human RAB31 / Rab22B Protein (His Tag) SDS-PAGE