

General Information

Gene Name Synonym

Adipocyte adhesion molecule; Coxsackie- and adenovirus receptor-like membrane protein

Protein Construction

A DNA sequence encoding the rat ASAM (Q8K1G0) (Met1-Met232) was expressed, fused with a polyhistidine tag at the C-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 18

Molecular Mass

The recombinant rat ASAM comprises 226 amino acids and predicts a molecular mass of 25.6 kDa. The apparent molecular mass of the recombinant protein is approximately 37 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat ASAM / CLMP Protein (His Tag) SDS-PAGE