

Human HSPD1/HSP60 (His & GST Tag) recombinant protein



Catalog Number: 502325

General Information

Protein Construction

A DNA sequence encoding the human HSP60 (NP_955472.1) (Leu 2-Phe 573) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human HSP60/GST chimera consists of 809 amino acids and has a calculated molecular mass of 88.7 KDa. It migrates as an

approximately 52-65 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human HSPD1 / HSP60 Protein (His & GST Tag)
SDS-PAGE