

# Human KEAP1 (His & GST & AVI Tag) recombinant protein



Catalog Number: 503395

## General Information

### Gene Name Synonym

Cytosolic inhibitor of Nrf2; INrf2; Kelch-like protein 19

### Protein Construction

A DNA sequence encoding the human KEAP1 (Q14145) (Gln 2-Cys 624) was fused with the N-terminal polyhistidine-tagged GST tag at the N-terminus linked by the AVI tag.

### Organism

Human

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Met

### Molecular Mass

The recombinant human KEAP1/GST chimera consists of 875 amino acids and has a calculated molecular mass of 99.2 kDa. It migrates as an approximately 85 kDa band in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile 20mM Tris, 500mM NaCl, 10% gly, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human KEAP1 / INRF2 Protein (His & GST & AVI Tag) SDS-PAGE