

Human CASP10/MCH4 (His tag) recombinant protein



Catalog Number: 518961

General Information

Gene Name Synonym

Caspase-10, CASP-10, Apoptotic Protease Mch-4, FAS-Associated Death Domain Protein, Interleukin-1B-Converting Enzyme 2, FLICE2, ICE-Like Apoptotic Protease 4, CASP10, MCH4

Protein Construction

A DNA sequence encoding the Human CASP10/MCH4 Val220-Ile480 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 25mM HEPES, 10mM DTT, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human CASP10/MCH4 (His tag) recombinant protein