

Human B3GAT3 (His tag) recombinant protein



Catalog Number: 519176

General Information

Gene Name Synonym

B3GAT3, Galactosylgalactosylxylosylprotein 3-beta-glucuronosyltransferase 3, Beta-1,3-glucuronosyltransferase 3, Glucuronosyltransferase I, GlcAT-I, GlcUAT-I, Gal beta-1,3-Gal-R glucuronosyltransferase,,, REF: C1005

Protein Construction

A DNA sequence encoding the Human B3GAT3 Glu72-Val335 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 150mM NaC, 2mM EDTA, 20% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human B3GAT3 (His tag) recombinant protein