

# Human MGMT (His tag) recombinant protein



Catalog Number: 519274

## General Information

### Gene Name Synonym

Methylated-DNA--protein-cysteine methyltransferase, 6-O-methylguanine-DNA methyltransferase, O-6-methylguanine-DNA-alkyltransferase, MGMT

### Protein Construction

A DNA sequence encoding the Human MGMT Met1-Asn207 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2 \mu\text{m}$  filtered solution of 20mM TrisHCl, 1mM EDTA, pH 8.0 .

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human MGMT (His tag) recombinant protein