

# Human HPGDS/GSTS/PGDS/PTGDS2 recombinant protein



Catalog Number: 519455

## General Information

### Gene Name Synonym

Hematopoietic Prostaglandin D Synthase, H-PGDS, GST Class-Sigma, Glutathione S-Transferase, Glutathione-Dependent PGD Synthase, Glutathione-Requiring Prostaglandin D Synthase, Prostaglandin-H2 D-Isomerase, HPGDS, GSTS, PGDS, PTGDS2

### Protein Construction

A DNA sequence encoding the Human HPGDS/GSTS/PGDS/PTGDS2 Met1-Leu199 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 200mM NaCl, pH 7.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human HPGDS/GSTS/PGDS/PTGDS2 recombinant protein