

# Bovine PDIA3 (His tag) recombinant protein



Catalog Number: 519626

## General Information

### Gene Name Synonym

PDIA3 protein, Protein disulfide-isomerase A3, 58 kDa glucose-regulated protein, Disulfide isomerase ER-60, ER protein 57

### Protein Construction

A DNA sequence encoding the Bovine PDIA3 Ser25-Leu505 is expressed with a 6His tag at the C-terminus.

### Organism

Bovine

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH7.5.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Bovine PDIA3 (His tag) recombinant protein