

# Human GSTP1/FAEES3/GST3 recombinant protein



Catalog Number: 519651

## General Information

### Gene Name Synonym

Glutathione S-transferase P;GSTP1;GST class-pi;GSTP1-1;FAEES3;GST3;

### Protein Construction

A DNA sequence encoding the Human GSTP1/FAEES3/GST3 Met1-Glu210 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM Tris, 150mM NaCl, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

## SDS-PAGE



Human GSTP1/FAEES3/GST3 recombinant protein