

Anti-S100A8 antibody



Catalog Number: 100010

Product name

Anti-S100A8 antibody

Immunogen

[Mouse S100A8 \(His Tag\) recombinant protein](#)

Specificity

Mouse S100A8

Antibody description

Rabbit polyclonal to S100A8

Preparation

Produced in rabbits immunized with purified, recombinant Mouse S100A8 (rh S100A8; P27005; Pro2-Glu89). S100A8 specific IgG was purified by Mouse S100A8 affinity chromatography.

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 μ g/mL

This antibody can be used at 0.1-0.2 μ g/mL with the appropriate secondary reagents to detect Mouse S100A8. The detection limit for Mouse S100A8 is < 0.039 ng/well.

IHC-P: 0.05-1 μ g/ml

Validations



S100A8 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry