

Anti-SFTPD / SP-D / SFTP4 antibody



Catalog Number: 100021

Product name

Anti-SFTPD / SP-D / SFTP4 antibody

Immunogen

[Rhesus SFTPD / SP-D / SFTP4 \(His Tag\) recombinant protein](#)

Specificity

Rhesus SFTPD / SP-D / SFTP4

No cross-reactivity in ELISA with Human SFTPD

Antibody description

Rabbit monoclonal to SFTPD / SP-D / SFTP4

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rhesus SFTPD / SP-D / SFTP4 (NP_001035283.1; Met1-Phe375).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions

IHC-P: 1-10 µg/mL

Validations



SFTPD / SP-D / SFTP4 Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of cynomolgus SFTPD in cynomolgus lung with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).