

# Anti-HMGB3/HMG4 antibody



Catalog Number: 100082

## Product name

Anti-HMGB3/HMG4 antibody

## Specificity

Human HMGB3

## Antibody description

Rabbit polyclonal to HMGB3/HMG4

## Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the human HMGB3, and purified by antigen affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.1-2 µg/mL

## Validations



HMGB3 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of human HMGB3 in human placenta with rabbit polyclonal antibody (0.5 µg/mL, formalin-fixed paraffin embedded sections).