

# Anti-P4HB antibody



Catalog Number: 100376

## Product name

Anti-P4HB antibody

## Immunogen

[Mouse P4HB \(His Tag\) recombinant protein](#)

## Specificity

Mouse P4HB

## Antibody description

Rabbit monoclonal to P4HB

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse P4HB (rM P4HB; NP\_035162.1; Met1-Lys506).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



P4HB Antibody, Rabbit MAb,  
Immunohistochemistry

Immunochemical staining of mouse P4HB in mouse brain with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded



sections).

P4HB Antibody, Rabbit MAb,  
Immunohistochemistry


Immunochemical staining of mouse P4HB in mouse liver with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded

# Anti-P4HB antibody



Catalog Number: 100376

---



P4HB Antibody, Rabbit MAb,  
Immunohistochemistry

Immunochemical staining of mouse P4HB in  
mouse intestine with rabbit monoclonal antibody  
(5 µg/mL, formalin-fixed paraffin embedded  
sections).

sections).