

Anti-B2M/beta-2 microglobulin antibody



Catalog Number: 100380

Product name

Rabbit IgG

Anti-B2M/beta-2 microglobulin antibody

Applications

IHC-P

Immunogen

[Rat B2M/beta-2 microglobulin \(His Tag\) recombinant protein](#)

Dilutions

IHC-P: 1-10 µg/mL


Specificity

Rat B2M / beta-2 microglobulin

Validations

Antibody description

Rabbit monoclonal to B2M/beta-2 microglobulin



B2M / beta-2 microglobulin Antibody, Rabbit MAb, Immunohistochemistry


Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Rat B2M / beta-2 microglobulin (rR B2M / beta-2 microglobulin; P07151; Met1-Met119).

Immunochemical staining of rat B2M in rat thymus with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).

Formulation

0.2 µm filtered solution in PBS



B2M / beta-2 microglobulin Antibody, Rabbit MAb, Immunohistochemistry

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Immunochemical staining of rat B2M in rat spleen with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).