

# Anti-S100A5 antibody



Catalog Number: 100383

---

## Product name

Anti-S100A5 antibody

## Immunogen

[Mouse S100A5 \(His Tag\) recombinant protein](#)

## Specificity

Mouse S100A5

## Antibody description

Rabbit monoclonal to S100A5

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse S100A5 (rM S100A5; P63084; Met1-Lys93).

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 0.01-2 µg/mL

## Validations



S100A5 Antibody, Rabbit MAb,  
Immunohistochemistry

Immunochemical staining of mouse S100A5 in mouse brain with rabbit monoclonal antibody (0.01 µg/mL, formalin-fixed paraffin embedded sections).